

**Alcohol
Wipes**



For the disinfection of general surfaces and medical devices



- Sani-Cloth 70% Wipes**
Disposable disinfection wipes for medical devices and general surfaces.
- Sani-Cloth Super Sani-Cloth Plus**
Disposable disinfection wipes for medical devices
- Sani-Cloth CHG 2%**
Disposable disinfection wipes for medical devices



Sani-Cloth 70% Wipes

Disposable disinfection wipes for medical devices.

Ingredients: Propan-2-ol 70 %v/v

- **200 Count Canister** Ref. **XP00159-VJT078** Packaging 6 case.
- **200 Count Refill Bag** Ref. **XP00198-VJT168** Packaging: 9 case.
Wipe size: 200X200mm Weight: 23g/m²
- **125 Count Canister** Ref. **XP00142-VJT057** Packaging 6 case.
- **125 Count Refill Bag** Ref. **XP00197-VJT101** Packaging: 9 case.
Wipe size: 130x185mm. Weight: 23g/m²
- **100 Individual Sterile Sachets Box** Ref. **XD24100-VJT061**
Packaging: 10 case. Wipe size: 200x128mm Weight: 37g/m²



Sani-Cloth 70% General Surfaces

Disposable disinfection wipes for general surfaces

Ingredients: Propan-2-ol 70 %v/v

- **200 Count Canister** Ref. **XP00285-VJT164** Packaging 6 case.
Wipe size: 200X200mm Weight: 23g/m²
- **200 Count Canister** Ref. **XP00284-VJT158** Packaging 6 case.
Wipe size: 220X200mm Weight: 23g/m²



Sani-Cloth Super Sani-Cloth Plus

Disposable disinfection wipes for medical devices

Ingredients: Propan-2-ol 70 %v/v, n-Alkyl dimethyl ethylbenzyl ammonium chloride 0.25%w/w, n-Alkyl dimethyl benzyl ammonium chloride 0.25%w/w.

- **200 Count Canister** Ref. **XP00141-VJT007**
Packaging: 6 case. Wipe size: 220X200mm Weight: 23g/m²
- **125 Count Canister** Ref. **XP00192** Packaging: 6 case.
- **125 Count Refill Bag** Ref. **XP00204** Packaging 9 case.
Wipe size: 130x185mm Weight: 23g/m²



Sani-Cloth CHG 2%

Disposable disinfection wipes for medical devices

Ingredients: Propan-2-ol 70 %v/v 2%w/v Chlorhexidine gluconate.

- **100 Individual Sachets Box** Ref. **XP00220-VJT103**
Packaging: 12 case. Wipe size: 150x162mm Weight: 25g/m²
- **200 Individual Sachets Box** Ref. **XP00314-VJT317**
Packaging: 10 case. Wipe size: 150x162mm Weight: 25g/m²



Features & Benefits

- Ready to use wipe.
- Deep well lid to store last dispensed wipe.
- 200ct slim canister for a better grip.
- Smear free.
- Brackets to fix on walls and/or vertical or horizontal poles (round or square poles).
- Colour coded red for all PDI alcohol based disinfectants & brackets to help identify formulation and provides compliance.
- Range of pack formats.

Compliance

- Pre-dosed wipe avoiding mixing error.
- Ready to use wipe therefore no activating period needed prior to use.
- Expressed liquid in compliance with its claims, ie. tested liquid.
- Fast contact time.
- 6 months PAO (applicable to canisters only).
- **FREE** Auditing along with product training available – Both supported by ATP system and tablets to manage collected data and ease access to e-training.
- Sani-Cloth CHG2% recommended for the disinfection of hubs and connection ports as per epic 2 and NICE guidelines – Provides residual disinfection.

Competitively priced

- Saves valuable time as ready to use.
- Refill bag is economical & environmentally friendly.



Alcohol Wipes

EFFICACY TEST	70% General CHG 2%	Super Plus
Bacteria EN 13727		
Pseudomonas aeruginosa	1min (c+d)	1min (c+d)
Staphylococcus aureus	1min (c+d)	1min (c+d)
Enterococcus hirae	1min (c+d)	1min (c+d)
DGHM surface test		
Pseudomonas aeruginosa	1min (c+d)	1min (c+d)
Staphylococcus aureus	1min (c+d)	1min (c+d)
Enterococcus hirae	1min (c+d)	1min (c+d)
Proteus mirabilis	1min (c+d)	1min (c+d)
Escherichia coli	1min (c+d)	1min (c+d)
EN1276		
Pseudomonas aeruginosa	30sec (c+d)	30sec (c+d)
Escherichia coli	30sec (c+d)	30sec (c+d)
Enterococcus hirae	30sec (c+d)	30sec (c+d)
Staphylococcus aureus	30sec (c+d)	30sec (c+d)
US EPA Guidelines (wipe test)-		
Pseudomonas aeruginosa ATCC 15442		5min (d)
Escherichia coli O157:H7 (ATCC 35150)		5min (d)
Staphylococcus aureus ATCC 6538		5min (d)
Klebsiella pneumoniae ATCC 4352		1min (d)
Salmonella		5min (d)
Methicillin resistant Staph aureus (MRSA)		5min (d)
Vancomycin Resistant Enterococcus faecalis		5min (d)
YEASTS EN 13624		
Candida albicans	1min (c+d)	1min (c+d)
(DGHM surface test)		
Candida albicans	1min (c+d)	1min (c+d)
FUNGI EN13624		
Candida albicans	2min (c+d)	2min (c+d)
ENVELOPED VIRUS DW/RKI		
Hepatitis B	15sec (c+d)	15sec (c+d)
Hepatitis C	15sec (c+d)	15sec (c+d)
Bovine viral diarrhoea virus (used as surrogate virus for Hepatitis C)	15sec (c+d)	15sec (c+d)
Vaccinia virus- Strain Elstree	15sec (c+d)	15sec (c+d)
HIV	15sec (c+d)	15sec (c+d)
US EPA Guidelines (wipe test)		
HIV - 1 Type 1		1min (d)
Influenza A Hong Kong		1min (d)
Vaccinia virus		1min (d)
NON ENVELOPED VIRUS US EPA Guidelines (wipe test)		
Adenovirus		1min (d)
Rhinovirus		1min (d)
TB- EN14348		
Mycobacterium terrae	1min (c+d)	1min (c+d)
DGHM surface test - Hygiene Nord		
Mycobacterium terrae	1min (c+d)	1min (c+d)
US EPA Guidelines (wipe test)-		
Mycobacterium bovis		1min (d)



*Subject to terms and Conditions

