



OXYGEN

FIRST RESPONSE


USER NAME		
ASSISTED BY		
DATE / TIME		
OXIMETER READING (SpO2)	INITIAL:	FINAL:
PULSE RATE (BPM)	INITIAL:	FINAL:
COMMENTS		



STEP 1
ATTACH TUBE TO
CANISTER AND
REMOVE SAFETY CAP



STEP 2
ATTACH TUBE TO MASK AND
PLACE OVER NOSE AND MOUTH.
PULL CORDS TO TIGHTEN



STEP 3
PRESS THE BUTTON ON THE
CANISTER TO RELEASE
OXYGEN AS REQUIRED